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Extraction parameters	
Cell type	Stainless steel
Cell size	10 mL
Oven temperature	60°C
Purge time	30 s
Nitrogen flow (gas assisted extraction)	10 mL/min per channel
Cell fill volume	50%
Solvent flow rate	1.0 mL/min
Extraction solvent	Acetonitrile:methanol (20:80 v/v%)
Extraction volume	≈22 mL
Pre-run rinse	10 mL, Acetonitrile:methanol (20:80 v/v%)
Extraction time (4 samples in parallel)	15 min

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Process 6 samples using the following steps:
1. Condition cartridge with 10.0 mL of 1% NH ₄ OH in methanol into solvent waste
2. Condition cartridge with 5.0 mL of 1% NH ₄ OH in methanol into solvent waste
3. Condition cartridge with 5.0 mL of 0.3M Formic acid into solvent waste
4. Pause and alert operator, resume when CONTINUE is pressed
5. Load 70.0 mL of sample onto cartridge
6. Pause and alert operator, resume when CONTINUE is pressed
7. Load 70.0 mL of sample onto cartridge
8. Pause and alert operator, resume when CONTINUE is pressed
9. Rinse sample bottle with 5 mL DI. Load 15.0 mL of sample onto cartridge
10. Pause and alert operator, resume when CONTINUE is pressed
11. Rinse sample bottle with 5 mL DI. Load 15.0 mL of sample onto cartridge
12. Pause and alert operator, resume when CONTINUE is pressed
13. Rinse sample bottle with 5 mL 1:1 Methanol: 0.1M Formic Acid. Load 15.0 mL of sample onto cartridge
14. Pause and alert operator, resume when CONTINUE is pressed
15. Dry Cartridge with gas for 0.8 minutes
16. Pause and alert operator, resume when CONTINUE is pressed
17. Manually rinse sample bottle with 5 mL 1% NH ₄ OH in methanol to collect
18. End

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