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Thermo Xcalibur 4.4采集

2.4.1 BioPharma Finder 4.1处理数据

The screenshot displays the BioPharma Finder 4.1 software interface, showing several configuration panels for peptide identification. The panels include:

- Basic Parameters:**
 - Enable HCP:
 - Protein Database: uniprot-bichi_jaomu_3A4922_and_USP2
 - Acquisition Type: High-High (high res MS1 and high res MS2)
 - Precursor Mass Tolerance: 20 ppm
 - Ions to Search: a ion, b ion, c ion, x ion, y ion, z ion, NL ion
 - Mass Range (MH+peptide mass range to analyze): 600.0000 to 5,000.0000
 - Identification Filters: E-Value Cutoff: 0.010000000 (all identification less than this value will pass filter)
- Protease Parameters:**
 - Select Protease: Trypsin
 - Protease Termini: Fully-digested, Semi-digested
 - Max Num Internal Miscalcavages: 0 0 1 2 (2 selected)
 - Enabled Decoy Search: Yes, No
- Advanced Search:**
 - Search by Full MS Only: Yes, No
 - Use MS/MS: Use All MS/MS
 - Maximum Peptide Mass: 11,000
 - Mass Accuracy (ppm): 10
 - Minimum Confidence: 0.80
 - Maximum Number of Modifications for a Decade: 1
- Disulfide Search:**
 - Perform Disulfide Bond Search: Yes, No
 - Allow Free Cys:
 - Maximum Number of Hits: 2,048
 - Maximum Number of Disulfide Bonds: 1
 - Maximum Number of Identical Chains in the Molecule: 2
- Variable Modification:**
 - Max # of variable modification per peptide: 1
 - 1. Variable Side Chain Modification: Deamidation (N) / +0.984 (N), N
 - 2. Variable Side Chain Modification: Deamidation (Q) / +0.984 (Q), Q
 - 3. Variable Side Chain Modification: Oxidation (MW) / +15.995 (M,W), M
 - 4. Variable Side Chain Modification: None

2.4.3 Spectronaut处理数据

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Protein Inference Workflow	From protein-db matching	Interference Correction	<input checked="" type="checkbox"/>
Digest Type	Specific	Only Identified Peptides	<input checked="" type="checkbox"/>
Enzymes/Cleavage Rules	<input checked="" type="checkbox"/> Trypsin/P <input type="checkbox"/> Trypsin <input type="checkbox"/> LysC <input type="checkbox"/> LysC/P <input type="checkbox"/> GluC <input type="checkbox"/> chymotrypsin(FLWY)	Exclude All Multi-Channel Interferences	<input checked="" type="checkbox"/>
Inference Algorithm	IDPicker	MS1 Min	<input type="text" value="2"/>
Toggle N-Terminal M	<input checked="" type="checkbox"/>	MS2 Min	<input type="text" value="3"/>
Exclude Duplicate Assays	True	Protein LFQ Method	Automatic
Generate Decoys	True	Proteotypicity Filter	None
Decoy Method	Mutated	Major (Protein) Grouping	by Protein Group Id
Preferred Fragment Source	NN Predicted Fragments	Minor (Peptide) Grouping	by Stripped Sequence
Decoy Limit Strategy	Dynamic	Major Group Quantity	Mean peptide quantity
Library Size Fraction	0.1	Major Group Top N	<input checked="" type="checkbox"/>
Machine Learning	Per Run	Max	<input type="text" value="3"/>
Precursor PEP Cutoff	0.2	Min	<input type="text" value="1"/>
Precursor Qvalue Cutoff	0.01	Minor Group Quantity	Mean precursor quantity
Protein Qvalue Cutoff (Experiment)	0.01	Minor Group Top N	<input checked="" type="checkbox"/>
Protein Qvalue Cutoff (Run)	0.05	Max	<input type="text" value="3"/>
Single Hit Definition	By Stripped Sequence	Min	<input type="text" value="1"/>
Exclude Single Hit Proteins	<input type="checkbox"/>	Quantity MS-Level	MS1
Pvalue Estimator	Kernel Density Estimator	Quantity Type	Area
		Data Filtering	Qvalue sparse
		Imputing Strategy	Global Imputing
		Cross Run Normalization	<input checked="" type="checkbox"/>
		Normalization Filter Type	None

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