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
Figure 2.Determining retention time windows

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Table 3. Initial demonstration of accuracy and precision

EPAMethod 1621	Prescribed limit	Experimental value
IDA (Section 9.2.1)	Recovery with 们± 20%	Recovery = 103.2%
IDP (Section 9.2.1)	RSD < 20 %	RSD = 4.47%

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Table 2. Comparison of relative standard error (RSE) and coefficient of determination ( $r^2$ )

	Linear			Quadratic			RSE limit
	No weighting	$1/A$	$1/A^2$	No weighting	$1/A$	$1/A^2$	
RSE*	38.76	20.47	7.25	13.22	5.94	5.75	20
$r^{2**}$	0.9967	0.9928	0.9900	0.9982	0.9977	0.9970	

$n = 7$   
\* RSE = Relative Standard Error, weighting = none,  $1/\text{area}$ , or  $1/(\text{area})^2$   
\*\*  $r^2$  = Coefficient of determination  
Calibration range = 2  $\mu\text{g/L}$  to 100  $\mu\text{g/L}$  based on a 100 mL sample


Figure 3. Calibration curve with RSE of 5.75

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
Figure 4. Method blank and 5 µg/L equivalent fluorine MDL standard

Table 4. Comparison of 3-day test of method blanks and 5 µg/L equiv. MDL standards

	Cindion C-IC and adsorption modules		Previous C-IC system reported in AN002748	
	Method blank F <sup>-</sup> (µg/L)	MDL standard* F <sup>-</sup> (µg/L)	Method blank F <sup>-</sup> (µg/L)	MDL standard* F <sup>-</sup> (µg/L)
Average over 3 days, X	2.17	4.57	1.23	4.76
Standard deviation (σ)	0.37	0.57	0.17	0.80
Calculated MDL**	—	1.7	—	2.5
Calculated ML(LOQ)	—	ML=5	—	ML=8

\* Average method blank subtracted  
 \*\*Section ii: MDL(s) = 3.14 \* σ<sup>40</sup>  
 Where MDL(s) is the MDL for spiked MDL standards, 3.14 is the Student's t-test value for n = 7, 99 percentile, and σ is the standard deviation.  
 Method limit (ML) = 3.18 \* MDL

Table 5. Sample results and recovery results of added PFHxS standard

Sample	Total fluoride ( µ g/L)		Recovered (%)
	Measured*	Added	
Municipal drinking water	1.81 ± 1.02	8.0	81.6 ± 1.6
		23.99	101 ± 16
Blended wastewater	4.86 ± 0.55	8.0	93.6 ± 14
		23.99	86.5 ± 5
Percolation pond water	0.29 ± 0.26	8.0	101 ± 5
		23.99	94.5 ± 7.3

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
Figure 5. Determination of PFAS in wastewater by AOF and C-IC

Table 6. Quality control checks for sample sequences

	Sample sequences	Total fluoride (µ g/L)		Recovered (%)
		Measured*	Added	
Method blanks	Drinking water	2.18 ± 0.28	--	--
	Wastewater	2.06 ± 0.09	--	--
	Pond water	2.59 ± 0.44	--	--
Low CV	Drinking water	4.45 ± 0.24	5.0	88.9
	Wastewater	5.05 ± 0.37	5.0	101
	Pond water	5.10 ± 0.12	5.0	102
OPR (n=1)*	Drinking water	23.4	25	93.6
	Wastewater	20.5	25	81.9
	Pond water	21.5	25	85.9
Medium CV	Drinking water	21.4 ± 1.15	25	85.5
	Wastewater	22.5 ± 1.02	25	90.1
	Pond water	25.20 ± 0.1	25	101
Method limits				80–120% OPR: 70–130%

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Term	Definition	Description
<b>AOF</b>	Adsorbable organic fluorine	Fluorocarbon and perfluorocarbon compounds in water or solvent that adhere to a sorbent like activated carbon or ion exchangers.
<b>CV</b>	Calibration verification standard	Check standards prepared from an inorganic fluoride standard that are from a different (second) source than the calibration standards. Added directly to combustion cups. Run at the beginning and the end of each sequence.
<b>GAC</b>	Granular activated carbon	Carbon material used for adsorbing organic compounds.
<b>IPR</b>	Initial precision recovery standard	Used to determine the expected process check value (OPR). Prepared by adding PFHxS standard to DI water. Processed through adsorption, combustion, absorption, and IC analysis.
<b>MB</b>	Method blanks	DI water processed through adsorption, combustion, absorption, and IC analysis. Indicates the lowest actual measurements possible and the baseline contamination of the C-IC and the filtration module. Always run at the beginning and the end of each sequence.
<b>MDL</b>	Method detection limit	The lowest reliable measurement, typically 3x S/N. In this test, the standard is processed through the adsorption on carbon, combusted, absorbed, and analyzed by IC.
<b>MDL(b)</b>	Method detection limit of blanks, n = 7	Defined as: (Mean of method blanks) plus (standard deviation x Student t-test factor).
<b>MDL(s)</b>	Method detection limit of 5 ng/mL PFHxS standard, n = 7	Defined as: (standard deviation) x (Student t-test factor).
<b>OPR</b>	Ongoing precision recovery standard	Check standard run at the beginning of each sequence. Prepared by adding PFHxS standard to 100 mL of DI water. Processed through the adsorption on carbon, combusted, absorbed, and analyzed by IC.
<b>PFHxS</b>	Perfluorinated sulfonate standard	Sodium perfluoro-1-hexanesulfonate standard. Used for MDL and QC check standards.
<b>RSE</b>	Relative standard error	Error from calibration line. Defines the best fit of a calibration curve.
<b>TSS</b>	Total suspended solids	Includes any loose particulates or precipitated solids in the water sample.

### Relative standard error formula

Definition of the relative standard error to determine instrument linearity

with

$$RSE = 100 \times \sqrt{\sum_{i=1}^n \frac{\left[ \frac{x'_i - x_i}{x_i} \right]^2}{n - p}}$$

$x_i$  = Nominal concentration (true value) of each calibration standard

$x'_i$  = Measured concentration of each calibration standard

$n$  = Number of standard levels in the curve

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