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
图1 Vanquish Analytical Purification Liquid Chromatography System

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Parameter		Value
Wash solvent		10% MeOH
Flush solvent		10% MeOH
Fraction collector compartment Temperature		5°C
Collection Time Frame		14.5-15.2min
		16.5-18min
		21-21.8min
Flush		Active
Max. tube volume		10ml
Minimum time for tube change		5s
Delay volume		900µl with delay capillary 1.0 x 1150 mm
Collection Path Mode		Horizontal
Collection Valve mode		Continue
Needle Positioning Mode		AboveVial
Use Safe Needle Height		Active
Wash Mode		Both
Wash Speed		100.0 µL/s
Wash Time		3.0 s
Rinse Mode		Both
Puncture Offset		0.0 mm
Peak Detection Options	Show advanced settings	
	Peak Start Threshold	570mAu
	Peak Start Slope	50 mAU/s
	Peak Start True Time	1.0s
	Peak End Threshold	570 mAU
	Peak End Slop	-50mAU/s
	Peak End True Time	1.0 s

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Parameter		Value	
Column		Hypersil Prep HS C18, 250 * 4.6 mm, 5 μm	
Eluent		A – 0.1%FA-H <sub>2</sub> O	
		B – MeOH	
Flow rate		1 mL/min	
Gradient	Time (min)	%A	%B
	0	92	8
	5	75	25
	13	70	30
	25	70	30
	28	0	100
	32	0	100
Autosampler temp.		5°C	
Column temp.		25°C	
Injection volume		100 μL	
Injection wash solvent		10% MeOH	
Detection		328nm、287nm、280nm	

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<sup>[1]</sup> Yujie Li, Qiangqiang Jia<sup>1</sup>, Min Zhang. et al Isolation of three glucaric acids from Leonurus japonicus Houtt. by using high-speed counter-current chromatography combined with semi-preparative high-performance liquid chromatography [J]. Journal of Separation Science. 2022;45:2140-2147;

<sup>[2]</sup> Dennis Koehler Principles of fraction collection using the Vanquish HPLC and UHPLC systems TN-72940

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热线 800 810 5118  
电话 400 650 5118  
[www.thermofisher.com](http://www.thermofisher.com)

**Thermo Fisher**  
S C I E N T I F I C