使用High-Resolution Accurate Mass Multi-Attribute Method (HR-MAM) 对NIST单抗的关键质量属性进行监控

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MS source setting	Value
Sheath gas	35 arb
Aux gas	10 arb
Sweep gas	0 arb
Spray voltage	3.5 kV
S-lens RF level	50
Aux gas temp.	250 °C
Capillary temp.	250 °C
Properties of Full MS - SIM	
General	
Runtime	0 to 72 min
Polarity	Positive
Full MS - SIM	
Resolution	140,000
AGC target	3e6
Maximum IT	200 ms
Scan range	300 to 1800 m/z
Properties of Full MS/dd-	
MS2 (Top5)	
General	
Runtime	0 to 72 min
Polarity	Positive
Default charge state	2
Inclusion	-
Exclusion	-
Tags	-
Full MS	
Resolution	140,000
AGC target	3e6
Maximum IT	200 ms
Scan range	300 to 1800 m/z
dd-MS2 / dd-SIM	
Resolution	17,500
AGC target	1e5
Maximum IT	250 ms
Loop count	5
TopN 5	TopN 5
Isolation window	1.2 m/z
Fixed first mass	-
(N)CE / stepped (N)CE nce	27
dd Settings	
Minimum AGC target	2.00e3
Intensity threshold	8.0e3
Apex trigger	-
Charge exclusion	Unassigned, 1, >8
Peptide match	Preferred
Exclude isotopes	On
Dynamic exclusion	8.0 s

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Database Parameters	
Protease	Trypsin (C-term KR)
Specificity	High
Fixed Modification	Carboxymethylation (C)
Variable Modifications	Deamidation (N)
	Deamidation (Q)
	Glycation (K)
	NH3 loss (NQ)
	Oxidation (MW)
	Lys (C-term)
	Gln-Pyro-Glu (N-term)
	N, O Glycans (CHO)

Table 3. Biopharma Finder 搜库条件.

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MS Detection	Value
Extracted Ion Chromatogram	<MS Default Detection Settings>
Detection Algorithm	ICIS
	Default Values
MS Settings	Value
MS No. of Decimal Points	5
Mass Tolerance	Manually Define Mass Tolerance, 8 ppm
Integration	Inhibit Integration for TIC Channel
MS Spectra Bunching	No baseline correction, 1 scan
Composite Scoring	Value
Pass Score if at Least	2 criteria passed
Fail Score if Less Than	1 criterion passed
Isotopic Dot Product	≥0.900
Mass Accuracy	≤5.00 ppm
Peak Apex Alignment	≤0.50 min

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Injection Name	MS Quantitation Lys Clipping(%)
190517_NIST_RA2_1_MS	89.148
190517_NIST_RA2_2_MS	89.150
190517_NIST_RA2_3_MS	89.150
190517_NIST_RA3_1_MS	88.782
190517_NIST_RA3_2_MS	88.332
190517_NIST_RA3_3_MS	89.006
Average	88.928
STDEV	0.325
CV	0.37%

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- Rogers, R., et al. Development of a quantitative mass spectrometry multi-attribute method for characterization, quality control testing and disposition of biologics. *MAbs*, 2015, 7, 881.
- Rogers, R., et al. A view on the importance of “multi-attribute method” for measuring purity of biopharmaceuticals and improving overall control strategy. *The AAPS Journal*, 2018, 20, 7.

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