

Evaluation of the AcquireX workflow for deciphering lipidome analysis of neutral lipids from whole insects using chromatography-based methods with high-resolution Orbitrap MSⁿ

Daniel Gachotte¹, Yelena Adelfinskaya¹, Jeffrey Gilbert¹, Jonathan Robeson², Adane Kassa², Reiko Kiyonami³, Yasuto Yokoi⁴ and David Peake³

Corteva, Indianapolis, IN, ²Corteva, Johnston, IA, ³Thermo Fisher Scientific, San Jose, CA, ⁴Mitsui Knowledge Industry, Tokyo, Japan



Agriculture Division of DowDuPont™



Introduction

Insect lipids are highly structurally diverse species that perform many important functions including storage of metabolic energy, contributing to the structure of membranes, protection against dehydration and pathogens, and circulating energy molecules. In addition, essential lipids such as sterols (precursors of hormones) and polyunsaturated fatty acids are only available from the diet. Understanding how lipid molecular species change in response to both diet and age is critical to define nutrient requirement and fitness relationship. Detailed studies of insect lipid composition however, have been technically challenging because of the complexity of lipid profiles. We compared the robustness of a faster LC high resolution LC-MSⁿ approach to perform non-targeted neutral lipid profiling experiments. Western Corn Rootworm larvae (*Diabrotica virgifera virgifera*) was chosen due to its economical impact estimated at 2 Billion USD cost associated with its control and corn production loss in the Americas. Three populations of larvae per condition were extracted and the lipids were mixed with commercially available labeled standards. Identification of neutral lipids present in the extracts of whole insects was achieved with two different chromatographic and acquisition methods.

Methods

Sample Preparation

Fresh frozen insect larvae were weighed and softened in hot IPA for 20 minutes to inhibit the lipases (Welti, Li et al. 2002). The samples were then homogenized using motorized glass rod and re-extracted for 10 minutes. After the addition of CHCl₃ and water the samples were agitated vigorously for 30 minutes at room temperature.

For tissue containing significant amount of neutral lipids a second extraction with CHCl₃:MeOH was performed under agitation at room temperature. The organic phase(s) was washed with 1 M KCl followed by water, filtered (PTFE disc 0.45 μm) and evaporated to dryness under N₂ prior to reconstitution in appropriate solvent containing Splash™ mix from Avanti Polar Lipids at 1/10 dilution. All samples were flushed with Argon and kept at -80°C.



Statistical Analysis

All the neutral lipids (Diacylglycerol, triacylglycerol, cholesterolesters) identifications resulting from ThermoFisher Scientific™ LipidSearch™ 4.2.9 software were determined in positive polarization modes. Only one adduct by lipid species with best quality score was used to compare four groups of samples. More than 12 measurements for each lipid was tabulated and used for visualization in charts. The tables were exported in Excel then JMP 12.2.0 (SAS Institute). In this poster only the species with highest degree of confidence (annotation Grade A or B) detected in 2/3 of the injections (biological replicates) were filtered to be compared using the t-test at p<0.01 confidence.


Overview Of The Lipidomic Workflow



LC: ThermoFisher Scientific™ Dionex Ultimate 3000™ RSLC system
Column: ThermoFisher Scientific™ Accucore™ C30 (2.1 x 150 mm, 2.7 μm, 45°C)
Condition 1: Mobile phase A: 60:40 AcN/H₂O + 10mM NH₄HCO₂, 0.1% HCO₂H
 Mobile phase B: 90:10 IPA/AcN + 10mM NH₄HCO₂, 0.1% HCO₂H
 Flow rate: 260 μL/min, Injection volume: 2 μL
Condition 2: Mobile phase A: AcN + 5mM NH₄OAc
 Mobile phase B: 75:25 IPA/MTBE + 5mM NH₄OAc
 Flow rate: 300 μL/min, Injection volume: 2 μL
Ref: Dionex HPLC-CAD lipid detection application note
Narveaz-Rivas, et al. J. Chrom. A 2016
MS: ThermoFisher Scientific™ Orbitrap Fusion MS with Optamax NG source
 AcquireX™ background exclusion mode described below:



Result: Sample Separation




Eluent profile of two LC separations compared in this study. The first condition could reduce acquisition time by half and eliminate the use of MTBE solvent for neutral lipid determination. Data was not collected during column re-equilibration time.

Splash mix was spiked in the samples at the concentration shown below. Chromatographic condition 2 was also tested with a Q Exactive Plus. Better reproducibility was obtained with the Fusion and both neutral lipid LC methods.


Class Components	Conc. (μg/ml)	Condition2			Condition1	
		RT (min)	QE_RSD (%)	Fusion_RSD (%)	RT (min)	RSD (%)
ChE 18:1(d7) Chol Ester	35	47.1	6.2	6.5	23.9	8
DAG 15:0-18:1(d7) DG	1	ND	NA	ND	16.9	5.2
TAG 15:0-18:1(d7)-15:0 TG	5.5	38.1	44	6.75	22.56	7.5
Ch Cholesterol (d7)	100	ND	NA	NA	NA	NA

The elution profile of neutral lipids is shown below. A more complete separation of DAG and TAG species is obtained using condition 2. The detection of DAG is better under condition 1.




Results: Diet and Age Effects


The sample set was divided into 4 classes with two developmental stages (1st and 3rd instar) and 2 diets. The samples were weighed to a target of 100 mg representing different number of larvae between stages (see pictures below). The Acquire-X method allowed Targeted Loss (fatty acid and ammonia) based identification at the MS³-level for TAGs and DAGs. A maximum of 415 TAG species were detected with chromatographic condition 2 and 60 DAG species with condition 1. Both chromatographic methods were necessary to collect a maximum of information. Peak areas were determined over ~4 orders of magnitude.



Age	Diet	
	Natural	Artificial
1st instar	Orange	Blue
3d instar	Green	Dark Blue



T-test results (p<0.01)			
Class	Lipid Species	AGE	DIET
TAG	415	98	295
DAG	60	23	29



Conclusions

Reproducibility measured by the variation of spiked standards representing 4 neutral lipid classes was significantly improved for two classes compared with the analysis performed on a Q Exactive. Chromatographic conditions did not influence the result. The Optamax® source and the Fusion instrument reduced variability.

The AcquireX strategy detected more than 358 and 432 neutral lipid species with chromatographic conditions 1 and 2, respectively. Both separations are complementary for DAG and TAG species. Diet conditions influenced the storage lipids levels to a greater extent than age difference.

The same sample set analyzed previously using the Q Exactive under the same separation conditions (Gachotte et al, poster 288550, ASMS 2017) gave 163 species annotated with similar confidence. The same acquisition with the Fusion™ instrument resulted in 415 species showing the power of this new instrument.

The AcquireX™ strategy allowed us to compile biological variability information (statistical analysis) and provided further annotation with MSⁿ (HCD and CID). This approach was specifically targeting DAG and TAG, but could accommodate other targeted lipid classes including other moieties after chemical derivatization.