Data Independent  
Acquisition

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**Thermo**  
SCIENTIFIC

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**Experiments**

- Full MS - SIM
- AIF
- Full MS / AIF
- Full MS / dd-MS<sup>2</sup> (TopN)
- Targeted-SIM
- Targeted-MS<sup>2</sup>
- Targeted-SIM / dd-MS<sup>2</sup>
- Full MS / AIF / NL dd-MS<sup>2</sup>
- DIA

**Properties of DIA**

**General**

- Runtime: 0 to 120 min
- Polarity: positive
- Default charge state: 2

**DIA**

- Resolution: 35,000
- AGC target: 2e5
- Maximum IT: auto
- Loop count: 30
- MSX count: 2
- MSX isochronous ITs: on
- Isolation window: 10.0 m/z
- Fixed first mass: -
- NCE / stepped NCE: 27

**Method editor — Inclusion List (modified)**

File	Edit	Help	Done	
Mass [m/z]	Formula [M]	Species	CS [z]	Polarity
1	835.62972			Positive
2	765.59788			Positive
3	915.66610			Positive
4	685.56150			Positive
5	515.48420			Positive
6	465.46146			Positive
7	705.57060			Positive
8	505.47965			Positive
9	955.68429			Positive
10	525.48874			Positive
11	815.62062			Positive
12	545.49784			Positive
13	985.69793			Positive
14	475.46601			Positive
15	745.58879			Positive
16	655.54786			Positive

**Method editor — Inclusion List (modified)**

File	Edit	Help	Done	
Mass [m/z]	Formula [M]	Species	CS [z]	Polarity
1965	525.48874			Positive
1966	405.43417			Positive
1967	975.69338			Positive
1968	635.53877			Positive
1969	825.62517			Positive
1970	675.55696			Positive
1971	765.59788			Positive
1972	465.46146			Positive
1973	875.64791			Positive
1974	955.68429			Positive
1975	735.58424			Positive
1976	425.44327			Positive
1977	575.51148			Positive
1978	805.61607			Positive
1979	905.66155			Positive
1980	815.62062			Positive

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## Large-Scale Proteome and Phosphoproteome Quantification by Data Independent Acquisition on Ultra-High Field Orbitrap Mass Spectrometers

Wei Zhang, Jing Li

Thermo Fisher Scientific, Shanghai, China

### Overview

1) Purpose: Data Independent Acquisition (DIA) has been applied for large-scale protein and PTM quantification on the Thermo Scientific™ Q Exactive HF and Orbitrap Fusion™ Tribrid™ instruments, to evaluate performance of the ultra-high field Orbitrap™ based DIA strategy for complicated sample analysis.

2) Methods: HeLa protein digest and enriched rat phosphorylated protein digest were analyzed on a Q Exactive HF MS and a Orbitrap Fusion MS by DDA and DIA methods, respectively. Data processing was performed on Skyline (version 3.1) using DDA database searching result (Proteome Discoverer version 1.4) as spectra library.

3) Results: A total of 21710 peptides corresponding to 3597 protein groups were extracted with high confidence (Q value  $\leq 0.01$ ) as well as good reproducibility from only 200 ng HeLa protein digest within 2-hour LC gradient. Moreover, A total of 9990 phosphorylated peptides were characterized and quantified with the same criterion. The results indicate the strengths of DIA on ultra high-field Orbitrap MS.

### Introduction

Orbitrap platforms based Data Independent Acquisition (DIA) is becoming increasingly popular, owing to its high sensitivity, selectivity and reproducibility compared to classic Data Dependent Acquisition (DDA) as well as other DIA techniques.

Meanwhile, the flexibility of Orbitrap MS permits the development of novel acquisition strategies as well as solving current puzzles in proteomics such as PTM quantification.

The latest ultra-high field Orbitrap technology combines the smaller size Orbitrap cell (20 mm ID) with 5 kV of the central electrode voltage, and realizes 1.8-fold scan rate increase at the same resolution or 1.8-fold resolution increase at the same frequency. State-of-the-Art performance and flexibility are thus accessible on ultra-high field Orbitrap based Q Exactive HF and Orbitrap Fusion mass spectrometers.

Herein, we utilized the ultra-high field Orbitrap MS for large-scale proteomic and phosphoproteomic DIA analysis, to evaluate performance and strengths of the ultra-high field Orbitrap based DIA strategy. The strategy shows outstanding quantification ability compare to current DIA techniques and provides a powerful tool for high-throughput proteomic sample application.

### Methods

#### 1) Materials

Commercial HeLa protein digest was purchased from Pierce™ (PN: 88328), with a loading amount of 200 ng for each run. Phosphoprotein sample was enriched from rat tissue and digested by trypsin, with a final loading amount of 700 ng for each run.

#### 2) Liquid Chromatography

Each sample was analyzed on a Thermo Scientific™ EASY-nLC™ 1000 nanoflow LC coupled to a Thermo Scientific™ Nanospray Flex™ ion source. The sample was loaded onto a nano-C18 column and separated at a flow rate of 300 nL/min with following gradient: 0–3 min, 3–7% buffer B (0.1% formic acid in acetonitrile); 3–95 min, 7–22% B; 95–113 min, 22–35% B; 113–116 min, 35–90% B; 116–120 min, 90% B.

#### 3) Mass Spectrometry

DDA performed 120K resolution MS scan and then triggered top 15 precursors (QE HF) or acquired under 3-second top-speed mode for 60K resolution MS/MS scans. The MS/MS AGC target value was set at  $1e6$  (QE HF) or  $5e4$  (OT Fusion) with 100–110 ms of max injection time (Figure 1A).

DIA collected high-resolution MS/MS spectra by sequential 25 amu step length (26 amu isolation window, 0.5 amu overlapping between windows) from  $m/z$  400 to 1200 for HCD fragmentation (28%). MS/MS ions were detected using 30K/60K of Orbitrap resolution,  $1e6$  (QE HF) or  $1e5$  (OT Fusion) of AGC target value and 85–100 ms of max injection time. MS scan was also performed before each DIA cycle (Figure 1B).

#### 4) Data Analysis

Database searching was performed on Thermo Scientific™ Proteome Discoverer™ software (version 1.4) with Uniprot human or rat protein database. Mass tolerance of precursor

and fragment were set at 10 ppm and 0.02 Da, respectively. Results were validated by Percolator with  $q < 0.01$  and then used as spectral library for DIA analysis.

DIA Data processing was performed on Skyline software (version 3.1). Three precursor (mono) isotopes and top six fragment ions in spectral library of each peptide were selected as transitions for peak extraction. Results were refined by mProphet using decoy sequences and then validated by Q value ( $< 0.01$ ).

### Results

#### 1) Data Independent Acquisition for Large-Scale HeLa Proteome Analysis

Three replicates were performed for each DDA and DIA methods, with a loading amount of 200 ng HeLa digest for each run. DDA and DIA were proceeded under identical sample and LC conditions to evaluate and compare results of the two strategy directly.

Figure 2 shows MS1 chromatograms of the DDA and DIA replicates. DIA limits cycle time at 3–3.5 seconds and thus maximizes sensitivity and accuracy by increasing AGC value, maximum injection time and resolution within the cycle time limitation. The zoomed-in chromatograms display scan point distributions of an example DDA and DIA peak. The cycle time of DIA is fixed and thus the chromatogram scan points are uniformly distributed, compared to that of DDA data.


FIGURE 1. Schematic comparison of data dependent (A) and data independent (B) acquisition methods on ultra-high field Orbitrap platforms.


FIGURE 2. Comparison of DDA and DIA raw data. (A) Comparison of MS1 base peak chromatograms. (B) Comparison of zoomed-in chromatograms. (C) Composite score distribution of targets and decoys of DIA data.

DDA ID results were imported into Skyline as spectral library for DIA data processing. mProphet was used to perform statistical analysis of targets and decoys, and generated Q value for result refinement (Figure 2C). DDA IDs with  $q < 0.01$  (Percolator) and DIA IDs with Q value  $< 0.01$  were extracted as high confident results.

A total of 23072 peptides corresponding to 3794 proteins were identified from only 200 ng HeLa digest by DIA, showing good sensitivity of the ultra-high field Orbitrap. Moreover, 94.1% of the peptide IDs and 94.8% of the protein IDs were obtained from all three replicates of DIA, showing much higher reproducibility than DDA, as 57.0% of peptides and 72.8% of proteins identified from all replicates of DDA (Figure 3 & 4).

### Conclusion

- 1) DIA has been applied and evaluated for large-scale protein and PTM quantification on the latest ultra-high field Orbitrap platforms.
- 2) A total of 21710 peptides from 3597 proteins were identified using DIA from only 200 ng HeLa digest, showing good sensitivity and reproducibility for quantification.
- 3) A total of 11341 phospho-peptides were identified using DIA and also displayed good sensitivity and reproducibility. Phospho-site ambiguity can be well distinguished, indicating the performance of DIA for PTM analysis.


FIGURE 3. Sensitivity and reproducibility of peptide ID results of DDA and DIA. (A) More than 90% of library peptides (DDA result) can be identified from DIA with Q value of all three replicates  $< 0.01$ . The ID reproducibility of DIA (C) is much better than DDA (B).

The variations of transition/precursor peak areas were also calculated and compared to evaluate quantification ability of DDA and DIA. The CV of 95.2% DIA peptides are below or equal to 20%, while only 74.0% DDA peptides have a  $CV \leq 20\%$  (Figure 5).


FIGURE 5. Comparison of peptide number with CV of transition/precursor peak area  $< 20\%$  (red) or  $> 20\%$  (blue). The quantification performance of DIA is much higher than DDA.


FIGURE 4. Sensitivity and reproducibility of protein ID results of DDA and DIA. (A) More than 90% of library proteins (DDA result) can be identified from all three replicates of DIA. The ID reproducibility of DIA (C) is also better than DDA (B).

### Data Independent Acquisition for Large-Scale Phosphorylation Analysis

The DIA strategy was further applied for phosphorylated protein sample to evaluate the performance of DIA for PTM analysis. A total of 11341 phospho-peptides were extracted with Q value  $< 0.01$  in all replicates. Figure 6 shows the DDA and DIA chromatograms as well as the DIA composite score distribution. Phospho-site ambiguity can be well distinguished due to high quality library spectra, good chromatogram separation and high sensitive DIA acquisition (Figure 7). DIA shows high ID reproducibility (Figure 8) and low peak area CV (Figure 9) as well in PTM analysis.


FIGURE 6. Comparison of phosphorylation sample data. (A) Comparison of MS1 base peak chromatogram. (B) Composite score distribution of targets and decoys of phosphorylation DIA data.


FIGURE 7. Sensitivity and reproducibility of phospho-peptide ID results. (A) More than 90% of library phospho-peptides (DDA result) can be identified from DIA with Q value of all three replicates < 0.01. The ID reproducibility of DIA (C) is much better than DDA (B).


FIGURE 8. Example peptides with identical sequence and mass but different phosphorylation site. (A) Correct peak assignment was obtained due to high data quality.

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类别	题目	作者	来源
Paper	Multiplexed MS/MS for improved data-independent acquisition	Egertson J D, Kuehn A, Merrihew G E, et al	Nat. Methods, 2013, 10: 744-746
Paper	Multiplexed peptide analysis using data-independent acquisition and Skyline	Egertson J D, MacLean B, Johnson R, et al	Nat. Protoc., 2015, 10: 887-903
Paper	DIA-Umpire: comprehensive computational framework for data-independent acquisition proteomics	Tsou C C, Avtonomov D, Larsen B et al	Nat. Methods, 2015, 12: 258-264
Paper	Extending the limits of quantitative proteome profiling with data-independent acquisition and application to acetaminophen treated 3D liver microtissues	Roland B, Oliver M, Bernhardt T G, et al	Mol. Cell. Proteomics, 2015, M114.044305
Paper	Multiplexed quantification for data-independent acquisition	Catherine E M, Alexander S H, Jarred W R, et al	Anal. Chem., 2015, DOI: 10.1021/ac503593d
Paper	Hybrid Data Acquisition and Processing Strategies with Increased Throughput and Selectivity: pSMART Analysis for Global Qualitative and Quantitative Analysis	Prakash A, Peterman S, Ahmad S, et al	J. Proteome Res., 2014, 13: 5415-5430
Paper	Quantitative Analysis of Targeted Proteins in Complex Sample Using Novel Data Independent Acquisition	Zhang W, Kiyonami R, Jiang Z, et al	Chin. J. Anal. Chem., 2014, 42: 1750-1758
Paper	Progress in Mass Spectrometry Acquisition Approach for Quantitative Proteomics	Zhang W	Chin. J. Anal. Chem., 2014, 42: 1859-1868
Poster	Characterizing Qualitative and Quantitative Global Changes in the Aging Heart Using a Novel HybridDA Method	Vogelsang M S, Prakash A, Sarracino D, et al	PN013 (USHUPO 2014) www.thermofisher.com
Poster	Large Scale Targeted Protein Quantification Using HR/AM Selected Ion Monitoring with MS/MS Confirmation on the Orbitrap Fusion Tribrid MS	Kiyonami R, Senko M, Zabrouskov V, et al	PN099 (USHUPO 2014) www.thermofisher.com
Poster	Large Scale Targeted Protein Quantification Using WiSIM-DIA Workflow on a Orbitrap Fusion Tribrid Mass Spectrometer	Kiyonami R, Patel B, Senko M, et al	PN 64085 (ASMS 2014) www.thermofisher.com

类别	题目	作者	来源
Poster	Evaluation of Data-independent Acquisition (DIA) Approaches for Spiked Peptides in HeLa Digest on Q-OT-qIT Mass Spectrometer	Zhang W, Kiyonami R, Jiang Z, et al	PN 64122 (ASMS 2014) www.thermofisher.com
Poster	Real-Time Qualitative and Quantitative Global Proteomics Profiling Using a Hybrid Data Acquisition Scheme	Schroeder T, Prakash A, Peterman S, et al	PN 64149 (ASMS 2014) www.thermofisher.com
Poster	Characterizing Qualitative and Quantitative Global Changes in the Aging Heart Using pSMART, a Novel Acquisition Method	Vogelsang M S, Prakash A, Sarracino D, et al	PN 64144 (ASMS 2014) www.thermofisher.com
Application Note	Large-Scale Targeted Protein Quantification Using WiSIM-DIA on an Orbitrap Fusion Tribrid Mass Spectrometer	Kiyonami R, Senko M, Zabrouskov V, et al	AN 600 www.thermofisher.com
Poster	Large Scale Targeted Protein Quantification Using HR/AM Selected Ion Monitoring with MS/MS Confirmation on a Novel Hybrid, Q-OT-qIT Mass Spectrometer	Kiyonami R, Senko M, Zabrouskov V, et al	PN W490 (ASMS 2013) www.thermofisher.com
Poster	Data-Independent Acquisition (DIA) Analysis on the Q Exactive Mass Spectrometer	Xuan Y, Kiyonami R, Egertson J, et al	PN Th229 (ASMS 2013) www.thermofisher.com
Poster	Proteomic Analysis of Acetaminophen Toxicity in 3D Human Liver Microtissues	Escher C, Ossola R, Gunness P, et al	http://cdn.f1000.com/posters/docs/299630429

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