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A

Sample	Description	Declared nicotine concentration (mg/mL)	Nicotine concentration (mg/mL)
a	Flavorless	0	<0.046
b	Flavored (branded)	0	<0.046
c	Flavored (branded)	0	<0.046
d	Flavored (vanilla)	0	<0.046
e	Flavored (mint)	0	<0.046
f	Flavored (branded)	6	5.5
g	Flavorless	6	7.3
h	Flavored (lemon)	12	12.2
i	Flavored (strawberry)	0	<0.046
j	Flavored (lemon)	0	<0.046

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1. United States Department of Health and Human Services. How Tobacco Smoke Causes Disease: The Biology and Behavioral Basis for Smoking-Attributable Disease: a report of the Surgeon General. Centers for Disease, Centers for Disease Control and Prevention, Atlanta, GA, USA, 2010.
2. Goriounova, N.A.; Mansvelter, H.D. Short-and Long-Term Consequences of Nicotine Exposure during Adolescence for Prefrontal Cortex Neuronal Network Function. *Cold Spring Harbor Perspectives in Medicine*, v.2(12), 2012.
3. World Health Organization Framework Convention on Tobacco Control signed in Geneva on 21 May 2003. [Online] <http://apps.who.int/iris/bitstream/handle/10665/42811/9241591013.pdf;sequence=1> (accessed Nov 15, 2018).
4. The Tobacco Products Directive (2014/14/EU). [Online] https://ec.europa.eu/health/sites/health/files/tobacco/docs/dir_201440_en.pdf (accessed Nov 15, 2018).
5. E-cigarettes: regulations for consumer products. Medicines and Healthcare products Regulatory Agency. [Online] <https://www.gov.uk/guidance/e-cigarettes-regulations-for-consumer-products> (accessed Nov 15, 2018).
6. Flora, J. W.; Wilkinson, C. T.; Sink, K. M.; McKinney, D. L.; Miller, J. H. Nicotine-related impurities in e-cigarette cartridges and refill e-liquids. *Journal of Liquid Chromatography & Related Technologies*, 2016, 39(17-18), 821–829.
7. Liu, X.; Joza, P.; Rickert, B. Analysis of nicotine and nicotine-related compounds in electronic cigarette liquids and aerosols by liquid chromatography-tandem mass spectrometry. *Beiträge zur Tabakforschung International/Contributions to Tobacco Research* 2017, 27(7), 154–167.
8. Aszyk, J.; Kubica P., Kot-Wasik, A.; Namiéśnik, J.; Wasik, A. Comprehensive determination of flavouring additives and nicotine in e-cigarette refill solutions. Part I: liquid chromatography-tandem mass spectrometry analysis. *Journal of Chromatography A* 2017, 1519, 45–54.
9. Herrington, B. J. S.; Myers, C.; Rigdon, A. Analysis of Nicotine and Impurities in Electronic Cigarette Solutions and Vapor. Restek, State College, PA, USA, 2015
10. CORESTA Recommended Method No. 62). Determination of Nicotine in Tobacco and Tobacco Products by Gas Chromatographic Analysis. 2005. [Online] <https://www.coresta.org/determination-nicotine-tobacco-and-tobacco-products-gas-chromatographic-analysis-29185.html> (accessed Dec 19, 2018).
11. Ghalop, V.V.; Kosmider, L.; Halquist, M.S. A Standardized Approach to Quantitative Analysis of Nicotine in e-Liquids Based on Peak Purity Criteria using High-Performance Liquid Chromatography. *Journal of Analytical Methods in Chemistry* 2018, Article ID 1720375.
12. Cooper, J.; Allen, C.; Cojocariu, C. Comprehensive chemical characterization of e-cigarette liquids using high resolution Orbitrap GC-MS. Thermo Scientific Application Note 10708, March 2019.
13. Davis, B.; Dang, M.; Kim, J.; Talbot P. Nicotine concentrations in electronic cigarette refill and do-it-yourself fluids. *Nicotine & Tobacco Research* 2015, 17(2), 134–141.
14. Trehy, M. L.; Ye, W.; Hadwiger, M. E. et al. Analysis of electronic cigarette cartridges, refill solutions, and smoke for nicotine and nicotine related impurities. *Journal of Liquid Chromatography & Related Technologies* 2011, 34(14), 1442–1458.

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