

20-July-2023

Dear Valued Customer,

Reference: Computers Specifications for Thermo Scientific GC and GC-MS Systems

A standard computer is included with all ISQ 7610 and TSQ 9610 series GC-MS systems. This is optionally available for TRACE 1600 series GC systems.

The standard computer is a Dell XE4 PC with Windows 10 Enterprise. Detailed PC specifications are as follows:

- Dell OptiPlex XE4 Minitower
- Win10 Enterprise Long-Term Servicing Channel (LTSC) IoT 2021
- Intel Core i7-12700 (12 core, 25MB, 20T, 2.1Ghz to 4.9Ghz, 65W)
- M.2 2230 PCIe NVMe 1TB Solid State Drive
- 32GB (4 x 8GB) DDR5 4800MHz
- Dell UltraSharp 24 Monitor U2422H
- Intel Integrated Graphics
- 8x DVD+/-RW Slimline optical disk drive
- US English (QWERTY) Dell KB-522 Wired Business Multimedia USB Keyboard
- Dell MS116 USB optical mouse
- Internal speaker
- 500W internal power supply unit (PSU), 92% efficient
- Intel vPro Enterprise
- Intel 2.5GbE i225 full height PCIe network adapter card (NIC)

Thermo Fisher Scientific greatly values our customers and if you have any questions on these products, please contact your local Thermo Fisher Scientific sales representative.

Best Regards,



Tim Anderson
Product Manager Single and Triple Quadrupole GC-MS
Thermo Fisher Scientific